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MEDICALTEST

DIAGNOSTIC TEST SPECIALIST

ENGLISH OVULATION TEST – INSTRUCTIONS

Ovulation test strips are samplers which contain a reactive for determining the presence of LH (luteinizing hormone) in urine samples. Ovulation tests help you maximise the chance of getting pregnant by determining the most fertile days of each cycle.

PROCEDURE

The ideal moment to perform the test is during the afternoon instead of the early morning. Perform the test preferably between 10 AM and 8 PM. In order to begin the test, open the sealed bag by tearing it along the groove. Take the LH test off the bag.

Carefully immerse the ovulation test vertically (with the arrow pointing down) in the urine container for 5 seconds, and afterwards place the test over a clean, dry, non-absorbent surface. Do not allow the urine level to exceed the line next to the lower part of the reactive test strip, otherwise the test will not be properly performed.

Wait for the coloured lines to appear. Depending on the LH concentration in the urine sample, positive results can be seen within just 40 seconds. However, up to 10 minutes are necessary in order to confirm negative results.


RESULTS

Within a period of three to five minutes, the two coloured lines will appear. Do not read the results after more than ten minutes have passed. In order to determine the result, compare the colour intensity; in other words, how light or dark the colour of the test line is in relation to the control line. In order to know whether a result is positive or negative, it is important to compare the colour intensity of both lines. This will tell us if an LH surge (which indicates ovulation) is taking place.

Positive increase in LH: if the test line has an equal or higher intensity (equal or darker) than the control line, this is a positive result and a good indication that the LH surge is taking place, and you are therefore ovulating.

Negative increase in LH: if the test line is a lower intensity (lighter) than the control line, or nothing appears on the test line, this means that the LH level of the sample is at or near normal levels, and an LH surge is not taking place. You are not ovulating.

Invalid result: if nothing appears on the control line within ten minutes, the result of the test is not valid and should be ignored. The control line does not appear if there is an insufficient amount of sample, or if the appropriate procedures have not been followed when performing the test. Repeat the procedure with a new ovulation test.

HOW TO RECOGNISE AN LH SURGE

After each test, you should know whether the LH surge is taking place. Commonly, six to ten days of tests are necessary to detect the LH surge, even though it is possible for additional testing to be necessary. LH surges generally last one to three days. Ovulation is more likely to take place one and a half days after the first day of the LH surge. If your test is negative, it is probable you are not having an LH surge. Remember that a test line that has a lighter pink colour than the control line only indicates that there is a very low level of LH in the urine.

Some of the most common explanations for a negative result are as follows:

Use of the first morning urine. The first morning urine should not be used to detect an LH surge.

LH concentration is too low to be detected. Continue performing tests on the following days until the LH surge is reached.

Testing stopped before the LH surge took place, and it should have been continued for a few additional days.

Testing is performed too soon or too late into the menstrual cycle (please, recheck the cycle calendar).

WARNINGS AND LIMITATIONS FOR THE TESTS

The test is only valid for in-vitro (external) diagnosis

Store the ovulation test in temperatures under 30 ° C. Do not freeze.

Do not use the ovulation test after the expiration date on the pack.

The ovulation test can only be used once.

Do not open the aluminium bag until you are ready to perform the test.

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